



The fully automatic Activity 885 Mark 2 is the successor of the popular Activity 885. The Activity 885 Mark 2 is up to 50 percent faster than its predecessor. The scanner convinces with its comfortable equipment that contains the software module secondDie next to a 36-month warranty. If secondDie is in use, the scanning process will be up to 40 percent faster. The Activity 885 Mark 2 achieves the data creation of a segmented 3-unit bridge in just 40 seconds. The reason for these amazing results is the speed update of the scanner software.

Of course the Activity 885 Mark 2 can capture articulated jaw models in perfect relation to the condyles. A fixator can be ordered additionally. Please take a look at our list of scanner accessories.

The scan data of articulated jaw models will be transferred to the virtual CAD through the patented system of scan fixators.

Your benefits at a glance:

- up to 50% faster than the previous model
- very high efficiency
- well-known smart optics quality „made in Germany“
- no licence fees

**Dimension ( B x H x T ) mm**

490 x 430 x 440

**Measurement field (X x Y x Z) mm**

80 x 60 x 85

**Accuracy (ISO 12836)**

6 µm

**Measurement Time Complete Jaw**

Scanning: 16 Sec.  
Postprocessing: 13 Sec.  
Total: 29 Sec.

**Measurement Time Single Tooth**

Scanning: 27 Sec.  
Postprocessing: 9 Sec.  
Total: 36 Sec.

**Measurement Time 3-unit bridge**

Scanning: 38 Sec.  
Postprocessing: 22 Sec.  
Total: 60 Sec.

**Weight**

32 kg

**Windows 7/8/10 64-Bit**

ScanFixator available

multiCase available

secondDie available

**System Tool Drawer****Power Supply Voltage**

100 – 240 VAC, 50/60 Hz

## Interfaces

USB

## Housing color

Anthracite  
RAL7016

## Compatibility with exocad® DentalCAD

### System requirements

#### Minimum

Windows 7 64-BIT  
Quadcore CPU, i3  
8 GB RAM  
USB 2.0 Port  
Approx. 40 – 100 GB free hard disk space  
Graphics card 1,024 MB RAM

#### Recommended

Windows 10 64-BIT  
Quadcore CPU, i7  
16 GB RAM  
USB 2.0/3.0 Port  
Approx. 80 – 150 GB free hard disk space  
Graphics card 2.056 MB RAM

### Deliveries include

1 3D scanner, 3 spacer plates, 1 locking plate, 1 hex key, 1 set of hex screws, 1 object holder, 1 mains cable, 2 USB cables, 1 measuring range template, 1 3D calibration model  
CD including scanner software, secondDie module, calibration data and operating manual (PDF)

The fully automatic Activity 885 Mark 2 is the successor of the popular Activity 885. The Activity 885 Mark 2 is up to 50 percent faster than its predecessor. The scanner convinces with its comfortable equipment that contains the software module secondDie next to a 36-month warranty. If secondDie is in use, the scanning process will be up to 40 percent faster. The Activity 885 Mark 2 achieves the data creation of a segmented 3-unit bridge in just 40 seconds. The reason for these amazing results is the speed update of the scanner software.

Of course the Activity 885 Mark 2 can capture articulated jaw models in perfect relation to the condyles. A fixator can be ordered additionally. Please take a look at our list of scanner accessories.

The scan data of articulated jaw models will be transferred to the virtual CAD through the patented system of scan fixators.

Your benefits at a glance:

- up to 50% faster than the previous model
- very high efficiency
- well-known smart optics quality „made in Germany“
- no licence fees